

Fast and reliable detection of radioactivity in food is achieved by the combination of a high performance handheld Gamma Personal Radiation Detector (PRD) and the portable Gamma Laboratory Kit.

Thermo Scientific Gross Gamma Food Monitor

Portable solution for food monitoring applications



- Portable, robust system in a transport case
- Minimal investment required
- Simple step by step operator's guide included
- RadEye PRD-S also to be used as a sensitive handheld gamma detector
- Battery power supply supports field operation
- Data logger for 1000 sample measurements
- PC interface via Infrared or optional Bluetooth
- Add-on scaler mode for regular RadEye PRD available (firmware upgrade)

The Gross Gamma Food Monitor consists of the RadEye PRD-S Gamma Pager and the specially tailored Gamma Laboratory Kit. Together they support the immediate requirements of the field assessment measurement program for contamination resulting from a nuclear accident.

The RadEye PRD-S gamma detector with its NaI(Tl) detector provides highly sensitive gamma measurements in a simple to use format.

The sample holder provides shielding from ambient radiation and an optimal counting geometry for the 180 ml samples.

Measurement results are provided as gross gamma rates in counts per second (cps) unit. Conversion factors translating the

'cps' results into Bq/l are provided in reference to the dominant elements (Iodine and Cesium) emitted by the Fukushima accident.

The system provides the capability to measure contamination levels for I-131 and Cs-134/137 to better than typical Derived Intervention Levels (DIL) in less than 5 minutes.

Additionally, the RadEye PRD-S can be used as an extremely sensitive handheld surveymeter. Personal radiation protection, checking of parcels for radioactivity and quick scanning of ravelers for I-131 incorporation (thyroid burden) are just some applications where the RadEye PRD-S delivers immediate and reliable measurement results.

Gross Gamma Food Monitor

To get an operational Gross Gamma Food Monitor, a Gamma Laboratory Kit has always to be combined either with a RadEye PRD-S or RadEye PRD-ER-S.

The sample holder of the Gamma Laboratory Kit contains 4 shielding modules that provide 4 pi shielding. Two sample containers, each holding 90 ml, are placed in a sandwich configuration around the RadEye PRD-S for optimized response. The optional Thermo Scientific 9 g Lutetium Test-Adapter provides the capability to perform QA checks on system operation.

Gamma Laboratory Kit (# 425069015)

Sensitivity


Typical Detection Limit (Measurement time: 300s)*	I-131: 140 Bq/l
	Cs-mix (50% Cs-134, 50% Cs-137): 115 Bq/l
Calibration Factor provided for	I-131, Cs-137, Cs-134
*Requires RadEye PRD-S or RadEye PRD-S-ER that is not included in the Gamma Laboratory Kit	

Carrying Case

Dimensions (approx.)	H 310 mm, W 360 mm, D 200 mm
Total weight (approx.)	20 kg



* Lu Test-Adapter positioning aid



3D assembly graphic



RadEye PRD-S (# 425067130) / RadEye PRD-ER-S (# 425067135)

Detector

Detector Type	Na(Tl) Detector with μ -Photomultiplier
Operation Modes	Ratemeter (cps, Sv/h, rem/h, R/h), Scaler (cps)
Measuring Range (Ratemeter mode)	RadEye PRD-S: 0.01 μ Sv/h - 250 μ Sv/h RadEye PRD-ER-S: 0.01 μ Sv/h - 100 mSv/h
Measuring Range (Scaler mode)	max. 65000 counts
Gamma Energy Range	60 keV - 1.3 MeV (dose rate), excellent detection from 30 keV
Firmware PRD-S / PRD-ER-S	Upgrade from regular RadEye PRD / PRD-ER available



© 2011 Thermo Fisher Scientific Inc. All rights reserved. Windows is a registered trademark of Microsoft Corporation in United States and other countries. Bluetooth is a trademark of Bluetooth SIG, Inc., Bellevue, Washington, U.S.A. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Results may vary under different operating conditions. Thermo Fisher Scientific makes no warranties, expressed or implied, in this product summary. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representatives for details. LIT110815_DB_Gross Gamma Food Monitor-e-V2.5

Europe, Africa, Middle East & Countries Not Listed

Fraunauracher Strasse 96 +49 (0) 9131 998-226
D 91056 Erlangen, Germany +49 (0) 9131 998-172 fax
customerservice.eid.erlangen@thermofisher.com

China

7th Floor, Tower West, Yonghe Plaza +86 10 8419 3588
No. 28 Andingm E. Street, Beijing, 100007 China +86 10 8419 3581 fax
info.eid.china@thermofisher.com

Singapore

11 Biopolis Way, Helios, Units #12-07/08 +65 6478 9728
Singapore 138667 +65 6478 9505 fax
info.eid.singapore@thermofisher.com

United Kingdom

Bath Road, Beenham, +44 (0) 118 971 5042
Reading RG7 5PR United Kingdom +44 (0) 118 971 2835 fax
customerservice.eid.beenham@thermofisher.com

USA, Canada, Mexico, Central & South America

27 Forge Parkway +1 (508) 553 1700
Franklin, MA 02038 USA +1 (800) 274 4212 US toll-free
info.eid@thermofisher.com +1 (508) 520 2815 fax

India

Plot No. C-327, T.T.C. Industrial Area, Pawne +91-22-41578800
Navi Mumbai 400 705, India +91-22-41578801 fax
info.eid.india@thermofisher.com

www.thermoscientific.com/rmp

Thermo
SCIENTIFIC